




QUICK *START* GUIDE

PRO 2.1





CARE & USE OF YOUR NEW TENNANT BIOMODULATOR® *PRO 2.1*

***PLEASE KEEP THIS GUIDE HANDY
AND REFER TO IT OFTEN***

Table of Contents

About this Guide.....	1
Batteries.....	2
Statement of User Responsibility.....	3
Contraindications.....	4
What is an Electrode.....	5
Care of Device.....	6
Can I do too much?.....	7
Powering on the Device.....	8
Accessory Port.....	9
Diagram – Front.....	10
Diagram – Back.....	11
Mode Selection.....	12
Left & Right Arrows.....	13
Modes.....	14-17
Assess Mode.....	18
Understanding Voltage Readings.....	19
Automatic Mode.....	20
How to Use Power.....	21
Fine Tune Power Setting.....	22
Graft of Power Curve.....	23

Table of Contents

Reaction Technology.....	24
Coefficients of Ongoing Reaction with Biofeedback.....	25
Step by Step Using Your Device.....	26-29
Treating Pain.....	30
Measure and Treat Each BioTerminal Point.....	31
Chart to Measure BioTerminal Points.....	32
Measuring Voltage Tips.....	33-34
Balancing Autonomic Nervous System.....	35-38
Activating Primary Power Circuit.....	39
Optional Attachments.....	40
Electrodes Attached with Wires.....	41
Improve Total Body Voltage.....	42-43
Chart Showing Where Black and Red Leads Go When Using Electrode Pads.....	44
Troubleshooting Voltage Readings.....	45
Simple Daily Treatment.....	46
Warranty.....	47
FDA Classification.....	48

About this Users Guide

- Please read all information in this guide and refer to it often. This material is intended to provide you with instructions to get you started using your device immediately. Senergy's daily protocol is very easy to follow and is the best way for all people to get started using their device.
- Please carefully read the **Contraindications** before using the device on anyone.
- You are encouraged to attend our health conference, as a way to confirm you have all the answers for using your device optimally. You will also gain a much greater understanding of how and why this device helps the body improve at the cellular level.
- You should **NOT** use your device in a professional therapy setting without having an in-depth understanding of proper use and contraindications.
- For assistance with your device, please reach out to the Senergy office at 972-580-0545 or visit our website, Senergy.us. Please note that our support does not include medical advice.

Batteries

- The Tennant BioModulator® operates using 2 AA alkaline batteries, which are included with the device. Please ensure you use new batteries of the same brand for optimal performance. **Avoid mixing old and new batteries, as this can lead to leakage. Battery corrosion can cause damage to your device.**
- When the ON/OFF switch is in the ON position, the device remains powered on, consuming battery life.
- To help conserve battery life, the device is designed to enter sleep mode after one hour. You can easily wake it up by pressing the (+) button.
- The device will alert you when the batteries are running low with a series of beeps or by reverting to Assess mode, displaying a low battery message. Please replace the batteries when you see this notification.
- If you have replaced the batteries and the beeping continues, please purchase a new pack of batteries and replace both with new ones.
- **It is very important to remove the batteries if you do not use your device regularly. Leaving batteries in the device for extended periods without use can lead to leakage, potentially damaging your device.**

Statement of User Responsibility

- You have received your **Tennant BioModulator® PRO 2.1** either as a licensed healthcare professional or with a prescription on file with Senergy Medical Group.
- Taking care of and using the device is entirely your responsibility.
- The examples shared by other users are meant to inform you and should not replace the guidance of your own physician or any other medical professional.
- These examples are intended to provide you with additional information to help you make informed decisions and **empower you to take control of your health.**
- **The power to manage your health lies within you.**



Contraindications

- These guidelines apply specifically to the use of the Tennant BioModulator® with the electrode placed directly on the skin. When using the Tennant BioTransducer®, these considerations do not apply.
- Please refrain from using it with a **pacemaker** unless you have the cardiologist's approval.
- It's important to avoid direct contact with any implanted electronics in the body.
- Please do not use on a **stroke patient** until any potential blood clot has had sufficient time to resolve.
- Refrain from rubbing the device over any vessels that may contain a **blood clot**, such as **tender varicose veins or the carotid arteries**.
- Please avoid using it when someone is under the **influence of alcohol or recreational drugs**.
- Do not apply the electrode directly to the eyeball.
- For liability reasons, please **do not use during pregnancy**.
- Please avoid using it on individuals who are very **sensitive to electromagnetic frequencies**, such as those who cannot be near a computer.
- **Do not use on a patient who expresses a desire not to have it used.**

What is an Electrode

- An electrode is a component that transmits current, allowing frequencies from the Tennant BioModulator® to enter the body.
- The metal on the back of the device serves as an electrode.
- All conductive attachments that connect to the device, such as foot plates, hand grips, Y electrodes, and electrode pads that adhere to the skin, are considered electrodes.
- Electrodes are conductive, which means they require skin contact for use, and you will experience a sensation when this contact occurs.
- If you are not feeling any sensation from an electrode attachment plugged into the side port of the device, it may indicate a faulty wire.
- The **Tennant BioTransducer®** is the **only attachment** that connects to the Tennant BioModulator® **that is NOT an electrode, so it does not require skin contact and you will not feel any current.**

Care of Device

Please keep in mind that your Tennant BioModulator®PRO 2.1 contains software and advanced circuitry similar to a computer, so it's important to handle your device with care..

- Please avoid turning it on if the temperature is below freezing.
- Kindly clean the electrode using only 70% isopropyl alcohol.
- We recommend not leaving your device in a hot car during the summer months.
- **Please avoid opening your device, as this may void your warranty.**
- Be cautious, as dropping the device on a hard surface can cause significant damage. We recommend using a protective carrying case.
- Unfortunately, if the device is dropped in water or if liquids are spilled inside, it will be irreparably damaged.
- The battery compartment is designed to slide on and off; please do not force the battery cover.
- Handle attachments gently when plugging them into the port, as excessive movement can lead to internal damage.
- Always unplug wires from the accessory port and store your device in a protective case for safekeeping.



Electrode corrosion can occur
If not cleaned regularly.



Can I do too much?

- You can use the device for extended therapy sessions whenever you feel the need. There may be moments when you're experiencing pain from a serious injury, and using it longer can support your healing process.
- It's essential to pay attention to your body's signals.
- If the current from the electrode feels too strong, please lower it. It should provide a gentle sensation, just a slight vibration.

Turn Tennant BioModulator[®] PRO 2.1 On

The [On/Off switch](#) is located on the side of your device. To power it on, please slide the switch up, and to turn it off, slide it down.



Accessory Port

Accessory port located in the top right corner of the unit.

Please use this port exclusively for authorized accessories.

Forcing anything into this port could void the warranty. so please handle it with care.

Remember to unplug the wire when it's not in use.



Diagram – Front



Diagram – Back



1. Built-in electrodes : Stimulation is applied through electrodes on the back of unit unless an accessory is attached.
2. Battery connectors: Only uses 2 alkaline 1.5V AA batteries.
3. Serial number
4. Calibration/programming connector – Do not damage or insert anything.

Very important to keep the battery compartment clean and dry.

Mode Selection with Tennant BioModulator® PRO 2.1

Menu Navigation
Keys- Use the
directional keys to
cycle through modes
and settings.

By pressing the down
arrow button, you'll be
able to see all the
modes arranged in the
same order as they
appear on the program
description sheet.



Select Key- Press
the O key to select
the desired option.

Left and Right Arrows

Advanced Function Keys

These buttons take you to sub-menu options for **advanced users**.
Ignore these buttons if you are a beginner.

The left arrow takes you to the Reaction (R) option when you are using one of the following modes;

Asses, Ten-8™, BT frequencies, CH frequencies, individual frequencies, and AVA™ modes.



The right arrow takes you to the option of fine tuning your power settings.

This lets you to increase the power using baby steps.

Tennant BioModulator® PRO 2.1 modes

- The following pages contain a list of operating modes, each programmed with specific frequencies and delivery patterns. You may observe that some modes offer a wide range of frequencies, while others provide a single frequency. When utilizing a mode with a single frequency, the goal is to deliver a concentrated frequency for shorter therapy sessions, typically lasting 5-10 minutes. In contrast, modes with multiple frequency changes can support therapy for extended periods, if necessary.
- Modes that incorporate Reaction technology are marked with (R). When you engage one of these modes, a number will appear in the center of the display screen, representing the voltage reading known as the initial reaction.
- Any mode can be paired with any attachment. An attachment is a device connected to the side of the Tennant BioModulator® to enhance the therapy experience. Common attachments include foot plates, hand grips, Y electrodes, and the Tennant BioTransducer®.

Tennant BioModulator® Pro - Modes


Mode Name / Frequency	Features – Types of things the program can help with
Assess (R) 60 Hz	Measurement – used to measure the BioTerminal points and other areas of the body, gives a (voltage reaction) reading. This is used with the electrode on the back of the device, placed firmly on skin.
Ten-8™ (R) 24- 30Hz	Pain & Injuries of Muscle, Tissue, Bone and Inflammation. - Used with BioTransducer® or any electrodes
Infinity™ 29 -121 Hz	Organs and Overall General Health - *Use this when you don't know what program to use.
Automatic (R) 29 - 121 Hz	Access (Measurement) + 1 Minute Therapy (Infinity Mode) This is helpful when you are measuring yourself and want to see if you can improve low areas with a quick 1-minute therapy of Infinity. Use with electrode on back of device.
SCR 147 (R) 147 Hz	Scars - This program is used to work on scars where energy blockage can be found, for the purpose of good energy flow through scar tissue. Use with BioTransducer® and soreness oil
BioTerminal® Step (BT Step) 8 -35 Hz	This program will do all the BioTerminal points for 30 seconds each, and repeats. Use with footplates and hand grips Can be used to balance body and emotions.
BioTerminal® Crown (BT-CR) (R) 15.18 Hz	Used to offer more voltage to this BioTerminal Point- Bring it back in balance.
BioTerminal® Head (BT-HD) (R) 13.53 Hz	Used to offer more voltage to this BioTerminal Point- Bring it back in balance.
BioTerminal® Neck (BT-NK) (R) 12.05 Hz	Used to offer more voltage to this BioTerminal Point- Bring it back in balance.
BioTerminal® Chest (BT-CH) (R) 10.74 Hz	Used to offer more voltage to this BioTerminal Point- Bring it back in balance.
BioTerminal® Abdomen (BT-Ab) (R) 10.14 Hz	Used to offer more voltage to this BioTerminal Point- Bring it back in balance.
BioTerminal® Pelvic (BT-PV) (R) 9.01 Hz	Used to offer more voltage to this BioTerminal Point- Bring it back in balance.
BioTerminal® Base (BT-BA) (R) 8.02 Hz	Used to offer more voltage to this BioTerminal Point- Bring it back in balance.
Visual Acuity (VAstep) 4 -1555 Hz and Infinity	Used to help with eye therapies. This program is used with the Transducer placed at the side of temple inline with the eye. See <i>Transducer Instructions for specifics on eye treatments.</i> *This is also a great program to use for all areas of the body.
Tennant Pain Management (TEN PM) 0.5 - 304 Hz	Pain & Injuries of Muscle, Tissue, Bone and Inflammation. Used with BioTransducer® or any electrodes
Solfeggio Musical Notes 174, 285, 396, 417, 528, 639, 741, 852, 963 Hz	This program has the frequencies of musical notes, resonates well with all areas of the body. (is used plugged into footplates and hand grips and with the Transducer) Overall General Health - * This is one of the best programs to choose when you don't know what program to use.
Harmonics (Schumann Waves) 7.83 -33.8 Hz	Grounding, like walking barefoot. Use for correcting polarity and when you have been around a lot of EMF, computers, Wi-Fi - Used with BioTransducer® for Polarity correction but can also be used for grounding with electrodes.

Mode Name / Frequency		Features
Delta Brain Waves	0.5 - 4.0 Hz	Deepest meditation and deep restorative sleep, Dreamless sleep. Healing and Rejuvenation. Delta waves suspend external awareness (located at back of head) - Occipital Lobe - Use with the BioTransducer®
Theta Brain Waves	3 - 8 Hz	Trance or hypnotic state more susceptible to hypnosis, daydream, twilight state, dream state, REM sleep. Subconscious, creativity and intuition. Improved memory, inspiration and spiritual connection, Improved physical healing, Release of beneficial hormones related to health and longevity. Reduction of mental fatigue, anxiety, and stress. Overall bliss and contentment (located above the ear) Temporal lobe - Use with the BioTransducer®
Alpha Brain Waves	7 - 12 Hz	'Frequency Bridge' between our conscious thinking (Beta) and subconscious (Theta) mind. known to help calm you down and promote feelings of deeper relaxation. Positive Thinking, Aid overall mental coordination, calmness, alertness, mind/body integration and learning, deep studying (located at top of the head) – Parietal Lobe – BioTransducer®
Beta Brain Waves	12 - 31.2 Hz	Normal Wakeful State, productivity, concentration, and alertness improved logic, reasoning, and critical thinking (located at front of head-) Frontal Lobe - Use with the BioTransducer®
Gamma Brain Waves	32 - 43 Hz	High Performance, perception is heightened, peak mental state. Improve memory and processing of information from different brain areas. Pass information rapidly and quietly. - Cerebellum - Lower back of the head- <i>Research suggests that people with learning difficulties or impaired mental processing may benefit.</i> - Use with the BioTransducer®
Chakra Step (CH step)	126 - 272 Hz	Cycles between all Chakra frequencies every 30 seconds, and repeats. - Use with the BioTransducer®
Chakra Root (CH Rt) (R)	194 Hz	Treat with BioTransducer® from the back-Located at Base between the legs (Perineum)
Chakra Sacrum (CH Sac) (R)	210 Hz	Treat with BioTransducer® from the back-Located two inches below the navel and center of the pelvis (backside)
Chakra Solar Plexus (CH Solrp) (R)	126 Hz	Treat with BioTransducer® from the back-Located along your spine and above the navel, in the upper abdomen (backside)
Chakra Heart (CH Ht) (R)	136 Hz	Treat with BioTransducer® from the back-Located back center of chest (backside)
Chakra Throat (CH Thr) (R)	141 Hz	Treat with BioTransducer® from the back-Located back of Neck
Chakra 3rd Eye (CH 3rd) (R)	221 Hz	Treat with BioTransducer® from back of head opposite of Forehead
Chakra Crown (CH Crn) (R)	172 Hz	Treat with BioTransducer® –Top of Head
Chakra 8th (CH 8th)(R)	272 Hz	Treat with BioTransducer® from back -Located arms length above your head (above the crown)
Blue Relax(R)	82 Hz	Similar to Assess but provides higher frequency output. Measures the body and gives a reaction reading.
Modulate	139 - 147 Hz	Acute pain and inflammation.
Stimulate	28 - 111 Hz	Scenar (FM/Var) Pain and inflammation- modulates both power output and frequency to reduce accommodation.
Deep Stim	30 - 121 Hz	Scenar (SW/Var) Pain and inflammation- multiple pulses and modulates both power output and frequency.
Blue Stim	22, 163, 22 Hz	Scenar - Pain and inflammation- this mode outputs multiple pulses and modulates frequency.
Acute	121 Hz	Scenar Russian technology—New Injury

Mode Name / Frequency		Features
Power 1	57 - 122 Hz	Toning face (when using on the face keep the power level below 35 and only treat for 1 - 2 minutes on each area, on clean dry skin) – Most people use the Face Electrode for toning the face.
Power Plus	57 - 122 Hz	Toning of muscles (has modulation) Use electrode pads for toning large muscle groups.
Power Advanced	57 - 122 Hz	Toning Body (higher intensity of pulses) When toning the larger muscle groups like the abs the electrode pads work well.
Power Intense	57 - 122 Hz	Toning Abs Muscles, Arms, Thighs (higher intensity of pulses) Use electrode pads for toning large muscle groups.
PG 2500	2300 - 2500 Hz	Highest Hertz- Outer energy field. Use with the BioTransducer®.
HGH	1725, 645, 1342 Hz	Human Growth Hormone - Most people use this with the BioTransducer®.
Repetitive Stimulate (R-Stim)	2500 Hz	Used by Physical Therapists
Fast T (R)	500 Hz	Used by Physical Therapists- Fast Twitch Muscle Fibers
Slow T (R)	250 Hz	Used by Physical Therapists- Slow Twitch Muscle Fibers
Lung (LU) (R)	824 Hz	Use by Acupuncturist-to treat meridian and to target organ for treatment- Used with BioTransducer® or any electrodes
Pericardium (PC) (R)	530 Hz	Use by Acupuncturist-to treat meridian and to target organ for treatment- Used with BioTransducer® or any electrodes
Heart (HT) (R)	497 Hz	Use by Acupuncturist-to treat meridian and to target organ for treatment- Used with BioTransducer® or any electrodes
Small Intestine (SI) (R)	791 Hz	Use by Acupuncturist-to treat meridian and to target organ for treatment- Used with BioTransducer® or any electrodes
Triple Burner (TE) (R)	732 Hz	Chinese Medicine all the organs; Upper - Heart, Lungs, Pericardium, , Middle -Spleen, Stomach, Gallbladder and Lower -Liver, Kidneys, Urinary Bladder, Small Intestine, Large Intestine and Uterus - . Use BioTransducer® at TB04 is located 'on transverse crease of dorsum of wrist in depression lateral to tendon of extensor digitorum communis.
Large Intestine (LI) (R)	553 Hz	Use by Acupuncturist-to treat meridian and to target organ for treatment- Used with BioTransducer® or any electrodes
Spleen (SP) (R)	702 Hz	Use by Acupuncturist-to treat meridian and to target organ for treatment- Used with BioTransducer® or any electrodes
Liver (LR) (R)	442 Hz	Use by Acupuncturist-to treat meridian and to target organ for treatment- Used with BioTransducer® or any electrodes
Kidney (KI) (R)	608 Hz	Use by Acupuncturist-to treat meridian and to target organ for treatment- Used with BioTransducer® or any electrodes
Bladder (BL) (R)	667 Hz	Use by Acupuncturist-to treat meridian and to target organ for treatment- Used with BioTransducer® or any electrodes
Gallbladder (GB) (R)	583 Hz	Use by Acupuncturist-to treat meridian and to target organ for treatment- Used with BioTransducer® or any electrodes
Stomach (ST) (R)	471 Hz	Use by Acupuncturist-to treat meridian and to target organ for treatment- Used with BioTransducer® or any electrodes
Advanced Mode (AVA)	0.5 - 1565 Hz	Where you can program additional frequencies to create additional therapy modes.
Memorized Mode (VAZ 4, 3, 2, 1)		Where you save the personalized modes, you created with the AVA program.
Emergency (ER)	24 - 30 Hz	Use in an Emergency for snake and spider bites as taught in the conference. Double select middle button takes you instantly to maximum power. - This is used with the electrode on the back of the device, placed firmly on skin.

Assess Mode

- **Assess** – “voltmeter” Reads the voltage when placed on the skin.
This is used to measure specific points on the body. Place the metal contact of the device on the skin to be measured. The voltage measurement is the Initial Reaction.




Understanding Voltage Readings

Your Initial Reaction Number

If your initial reaction number is between

- 5 – 14 this is a low range and needs improvement
- 15 – 24 this range is still slightly below normal
- 25 – 34 this is the normal range
- 35 – 49 this is normal range for athletes
- 50 and higher indicates healing or some inflammation



Automatic mode

- This operates in **two phases**: (1) It first measures your voltage, and then (2) the Infinity™ therapy mode starts automatically and lasts for one minute. You will hear a series of beeps when the therapy session concludes. It will then take another reading, allowing you to see if there has been any improvement. You can either move to treat another area or stay in place for the next treatment cycle.
- **Use Automatic when you wish to measure the voltage and begin treatment right away.**
- The aim is to elevate your initial reaction to above 25, which should occur within 2 to 3 minutes.
- If there is no improvement in your initial reaction after 2 to 3 minutes, switch the device to treatment mode and continue for several minutes.
- **This device is not designed for use with electrode pads.**
- **It is not intended for treatment durations longer than a few minutes, as it resets to measure the body every minute.**

Setting Power Level

- Please press and hold the (+) button to gradually increase the power.
- Adjust the power to a level that feels gentle and comfortable for you. You should feel a slight tingle.
- It's always safer to start with a lower power setting than to use too much.

Using too much power may cause the cells to shut down.


(This is not harmful but it is counterproductive.)

- Allow children to take charge of their own power.
- In Ten-8™ mode, the sensation will remain steady in intensity.
- In Infinity™ mode, the experience starts gently and builds in intensity over time, moving from low to high levels. Please take care and increase the power gradually. We suggest beginning with the power level in Ten8 before switching to Infinity, allowing you to set your maximum power level from the start.

Fine-Tune Power Setting

- This option is designed for individuals who may find the normal power settings too intense. It gradually increases the power in small increments.
- To adjust your power setting, press the right arrow followed by the up arrow. Your options include cosmetic, soft, sensitive, and ultra soft. Please confirm your choice using the middle button.
- After making your selection, the power level will reset to 1.
- This feature helps ensure that someone is not startled by a sudden increase in power.
- The default power setting is normal. Each time you turn off your device, it will revert to the normal power setting when you turn it back on.
- The next page includes a graph to illustrate the fine-tuning of the power level.

Graph of Power Curve



Reaction Technology D and Z

- The device shows an initial reaction in modes marked with (R).
- This initial reaction refers to the voltage measurement discussed in the Assess mode, which is displayed at the top center of the screen.
- You will also notice a D and Z appearing in the bottom corner of the screen. The D signifies that a dynamic change is occurring. The Z indicates that there is no change.
- The D and Z are there solely to provide you with information, which can be useful in understanding your treatment's progress.
- If you see a Z first, it means there is no dynamic change, and your body may require more treatment.
- Many users tend to overlook the D and Z when treating pain or when following a protocol with a set duration.

Coefficients of Ongoing Reaction(X) and Initial Reaction (Y)

This function is for advanced users, and you will need to reference the user manual for additional information.

- Reaction mode displays; Ongoing Reaction, Coefficient of Form, D and Z progress.
- Reaction mode is available for Assess, Blue Relax, Ten-8, BT modes, CH modes, individual modes and AVA™ modes, as indicated with (R)
- To enable the Reaction option, use the < left arrow key.
- To get back out of Reaction use the > right arrow key.
- If $X=Y$, there is no dynamic change, then the body needs more treatment.
- If $X<Y$, indicates that the body may need more treatment.
- If $X>Y$, the body is approaching a normal response to the stimulation.

Selecting A Mode

Navigation Keys: You can use the up or down buttons to scroll through the modes. By pressing the down arrow button, you'll be able to see all the modes in the same order as they appear on the program description sheet.



Once you've reached your desired mode, please select it using the middle button.

Put The Device On Your Skin And Press The + Button Until You Feel A Slight Tingle

- Gently move the device across your skin to familiarize yourself with its sensation.
- It's essential to feel at ease while using your BioModulator.
- Smoothly glide it over your skin.
- Move it up and down your arms, legs, torso, hands, and feet.



Using your device in this manner offers a valuable benefit.

Do Not Be Overly Concerned About Making Mistakes

- Some users may initially worry about making mistakes. However, if you follow the instructions, there's no need for concern.
- The Tennant BioModulator ® is designed for general use and does not require precise application.
- Even if you apply the therapy to a different area than intended, it is safe and will often still provide benefits.
- There are many effective ways to deliver additional energy, **so please don't stress about making mistakes.**



Everyone begins their journey with the Daily Protocol.

- This protocol addresses every area of balancing and charging up the Body.
- Please make sure to follow **your daily protocol steps consistently for all conditions.** After that, use your Tennant BioTransducer® to concentrate on the areas that need support or are experiencing pain.

Treat Pain As Long As Needed



- There are various options available for managing pain, and Ten PM is an excellent choice. It offers frequencies ranging from 0.5 to 304 Hz and delivers different pulse patterns. You can truly appreciate the differences in modes when you apply them directly to your skin using an electrode. Examples of electrodes include foot plates, hand grips, the Y electrode, and adhesive electrode pads.
- The Tennant BioTransducer® is another effective tool for addressing your pain. Often, it can provide deeper relief, especially in cases like a broken bone or swollen areas.
- It's beneficial to explore multiple modes and application methods for pain relief.
- Ten-8 and Ten PM are popular modes, but there are many others, such as Modulate, Stimulate, Deep Stim, Blue Stim, and Infinity, all of which can assist in alleviating pain.
- Once you've selected a mode, the next step is to determine how you want to apply the therapy. You can connect your Tennant BioTransducer® and direct the light towards the area of pain.
- The therapy can be administered for as long as you need.

Measuring and Treating BioTerminal® Points

- Measuring the body helps us identify areas that may need support.
- When taking voltage readings, please place the electrode on your device at each point indicated in the image on the next page, ensuring firm contact with the skin.
- If the initial reaction is below 25, this suggests low voltage.
- Ensure the skin is clean and dry, free from makeup or moisturizers.
- Use the **Assess or Automatic modes** to take measurements; **your voltage measurement will be your initial reaction number.**
- To **measure and treat simultaneously, use Automatic mode and continue treatment until the initial reaction reaches 25.**
- **You can treat all Bio Terminals at once by utilizing the foot plates and hand grips to energize the entire body.**
- Daily measurements are not necessary. Only the first measurement is significant; if you measure the same spot repeatedly, the results may vary each time.

For educational purposes
only

Measure Voltage on each BioTerminal Point



You can make copies of this chart for future use

Avoid Measuring Voltage When

- The skin is hot and sweaty
- The skin is very dry
- You have not had at least 8oz. of water to drink in the last 30 min.

If you take readings under these conditions, they will not be accurate

- **If the skin is too dry the reading will be lower**
- **If the skin is sweaty the reading will be higher**

Only Measure Voltage Once A Day

- When measuring your points, please ensure the electrode is placed firmly on the spot.
- The first reading you obtain is the most significant.
- Keep in mind that your voltage may vary from day to day.
- During exercise, your voltage tends to be higher as your muscles function like rechargeable batteries, storing voltage.
- It's essential to drink at least 8oz. of water 30 minutes prior to measuring your voltage.
- Remember, taking measurements is meant to help you gather information; it's not necessary to do this every day.

Balancing the Autonomic Nervous System

It's important to correct imbalances in the body.

The effects of balancing the autonomic nervous system can be both visible and subtle, yet they are certainly felt. This treatment is crucial because an imbalanced autonomic nervous system impacts every aspect of your well-being. It plays a significant role in chronic pain and chronic illness.

Here are some physical effects you may notice:

- There may be noticeable asymmetry in the two sides of the face.
- The jaw might shift to one side, potentially leading to TMJ issues.
- Sinus and nasal obstructions can contribute to snoring.
- The C1 and C2 vertebrae may misalign, resulting in headaches, neck pain, and migraines.
- The spine may curve, which can cause disk extrusion.
- One shoulder might be elevated compared to the other, making one arm appear shorter.
- The pelvis may rotate, leading to lower back pain, disk extrusion, and possibly causing one leg to seem shorter than the other. This can increase stress on one hip, knee, and ankle, leading to wear in those joints.
- The cranio-sacral pump may become locked, causing the nervous system to rely on stagnant cerebrospinal fluid, which can reduce overall function; however, restoration can help correct these issues.
- This can also contribute to improvements in blood pressure and eye pressure.
- Chronic headaches may find relief.


Balancing the Autonomic Nervous System

This is done by using the Tennant BioModulator® in the Infinity mode.

1. Place the Y electrode on the side of the neck half-way between the shoulder and the ear.(See image on next page.) device should be in Ten-8™ mode to set the power at a comfortable level. Switch to the Infinity™ mode to do the therapy.
2. You will need to push firmly into this muscle, as you are trying to reach the nerve behind the muscle.
3. When you find the right spot you will know, because it feels different. You should feel a tingling sensation along the back of the ear. It may take time to find the right spot so be patient.
4. Hold unit in place for 3 minutes.
5. Repeat treatment on other side of neck.
6. Once both sides are complete, turn head left to right, approximately 10 times each side.
7. Rotate shoulders front and back, approximately 10 times each way.
8. Stand and rotate at the waist, approximately 10 times.
9. Walk around room or area for a short distance.

This Image Shows Where to Place the electrode for Balancing the Autonomic Nervous System

Fixing Your Bowling Ball



Balancing the Autonomic Nervous System

We want to take a moment to emphasize the importance of the bowling ball technique. This is a skill that requires practice, and it's completely normal for it to take time to master.

Balancing the autonomic nervous system is a daily commitment, and it's essential for maintaining that balance.

We've noticed that those who are just starting out often find it easier to use the Infinity mode rather than the automatic setting. This approach allows you the freedom to explore beyond just one minute, as it may take that long to discover the right spot. When you do find it, you'll recognize it by a distinct feeling. You should start to notice sensations radiating in different areas, including a tingling sensation along the back of the ear. Once you experience that unusual tingle, you'll know you're on the right track. Keep the BioModulator in place for at least three minutes after feeling that sensation; the longer you maintain it, the more relaxed you will become.

Please be cautious with the intensity in Infinity mode, as it modulates at a higher level. Mastering this technique is crucial for effectively managing daily stress, and we're here to support you in this journey.

Activating Primary Power Circuit



Place the Y electrode at the base of the skull. Treat using the Ten-8 mode for 3 minutes

Optional Attachments

The unit features an Accessory Port on the right side. There are various attachments available for the **Tennant BioModulator®** that can help simplify the treatment of specific areas of the body.



• **Scalp Electrode**

- *For use in hair on scalp, chest, etc.*
- *Can also be used on animals*



• **Face Electrode**

- *Necessary for taking readings of voltage on ones, own face or back.*
- *For general use in small areas*



• **Pencil Electrode**

- *For use on acupuncture points*



• **Y Electrode**

- *Provides an 11-inch extended reach.*
- *Useful for practitioners too*

Always unplug any attachment when not in use.

Electrode Pads Attached With Wires



The electrode pads are connected by a wire that splits into two parts: one end is red and the other is black, which plugs into the pads. Depending on the wire you choose, you can use either 2 or 4 pads. This is a convenient way to manage pain while you're on the move.



Use this wire to connect 2 pads




Use this wire to connect 4 pads at once

Improve Total Body Voltage Charge the Muscle Battery Packs

- This is done by using the foot plates and hand grips at the same time. This is recommended for people to do everyday to charge up the whole body. Please charge after you have done the balancing steps of the daily protocol.
- You would use the **Solfeggio** mode for about 30 minutes a day.
- The image on the next page is showing the foot plates and hand grips. These attachments are highly recommended for everyone to use to charge up their body.

Foot Plates and Hand Grips Charge up the Muscle Battery Pack



Place Black Lead Wherever You See “Black”

Place Red Lead Wherever You See “Red”

If you are using the foot plates the right foot uses the red plate and left is black and hand grips hold the red in the left hand and the black in the right hand.


- This image illustrates the optimal placement of electrode pads to ensure effective voltage flow when using them on various body quadrants, such as arms and legs.

- If you're experiencing pain, please place one pad on the affected area, connecting the red and black wires as shown in the image. Position the other pad on the corresponding site.

- The pads can be arranged from left to right or front to back.

- This aligns with the polarity in our bodies.

POLARITY CHART



Troubleshooting Voltage Readings

- If you can not get a reading; resolve these issues first
 - Person is *dehydrated* – *drink more water*
 - *Reading takes longer to display (may take up to 10 seconds) when voltage is very low*
 - *At times, the voltage can be too low to be able to get a reading because of a serious condition – start a treatment by using **Ten-8™** and **Infinity™** modes. You will notice when voltage is especially low you will need to increase the power, and as voltage is improved in the area, you will need less power, (you continue to adjust the power as needed to keep it at a comfortable level)*

Simple Daily Treatment

- Drink 8oz of pure alkaline water before doing your protocol.
- Do your daily protocol as outlined in separate handout.
- Balance your autonomic nervous system whenever you feel stressed.
- Now that you have balanced your body and provided a base of energy for your body you can work on any areas of pain or chronic illness.
- An easy way to work on chronic pain or illness is with the Tennant BioTransducer®, simply point the light at any area of concern.
- Please call Senergy at 972-580-0545 to order a Tennant BioTransducer® if you don't have one or visit Senergy.us

Two-Year Limited Manufacturer Warranty

- The warranty begins on the date of purchase and covers manufacturer defects only.
- The warranty becomes invalid upon any of the following:
 - More than two years has passed since the purchase date.
 - The device was damaged or broken due to abuse, misuse or improper use.
 - Dropping it, throwing it, leaving it in extreme heat or cold, etc.
 - Water or chemical damage.
 - Broken seal or case; trying to open device.
 - Shorting of electrodes.
 - Violent vibration.
 - NOTE: violent vibration may damage components and reduce the effectiveness of the device.
 - To receive an additional year of warranty at no cost, please refer to your owner's manual under the Warranty section.

FDA Classification

Important Notice

FDA Classification: 882.5890 Neurology transcutaneous electrical nerve stimulator for pain relief

Identification for Use: Symptomatic relief and management of chronic, intractable pain

Adjunctive treatment in the management of post-traumatic surgical and post-traumatic pain

Relaxation training and muscle reeducation

Caution:

