

Tennant BioModulator Pro - Modes

Program	Feature
Assess	Measurement
Infinity	Organs, Use this when you don't know what program to use
Ten-8	Muscle, Tissue, Bone, Injury, Pain
Automatic	Measurement + 1 Minute Therapy
SCR 147	Scars
BioTerminal Step (BT Step)	This frequency will do all the BioTerminal points in order over and over (used plugged into footplates and hand grips)
BioTerminal Crown (BT-CR)	Used to offer more voltage to this BioTerminal Point- Bring it back in balance
BioTerminal Head (BT-HD)	Used to offer more voltage to this BioTerminal Point- Bring it back in balance
BioTerminal® Neck (BT-NK) (R) *	Used to offer more voltage to this BioTerminal Point- Bring it back in balance
BioTerminal® Chest (BT-CH) (R) *	Used to offer more voltage to this BioTerminal Point- Bring it back in balance
BioTerminal® Abdomen(BT-Ab)(R) *	Used to offer more voltage to this BioTerminal Point- Bring it back in balance
BioTerminal® Pelvic (BT-PV) (R)	Used to offer more voltage to this BioTerminal Point- Bring it back in balance
BioTerminal® Base (BT-BA) (R) *	Used to offer more voltage to this BioTerminal Point- Bring it back in balance
Visual Acuity Step (VA step)*	Used to help with eye therapies
Tennant Pain Management (TEN PM) *	Used for Pain, injuries, inflammation
Solfeggio	Frequency is like musical notes, resonates well with body. Can hear different tones
7.83 Harmonics (Schumann Waves)	Grounding, like walking barefoot. Use when out of sorts

Program	Feature
Delta Brain Waves 0.5 - 4.0 Hz	Deepest meditation and dreamless sleep. Delta waves suspend external awareness (located at back of head) –Occipital Lobe
Theta Brain Waves 3-8 Hz	learning, memory, and intuition (located above the ear) Temporal lobe
Alpha Brain Waves 7-12 Hz	aid overall mental coordination, calmness, alertness, mind/body integration and learning, deep studying (located at top of the head) –Parietal Lobe
Beta Brain Waves 12-31.2 Hz	Normal Wakeful State, (located at front of head-) Frontal Lobe
Gamma Brain Waves 37-43 Hz	Improve memory (high frequency, like a flute), and relate to simultaneous processing of information from different brain areas. Pass information rapidly and quietly. (located at cerebellum)- Lower back of the head
Chakra Step (CH step) (R)	Used by Acupuncturists cycles between all BioTerminal points over and over
Chakra 8th (CH 8th)(R)	Used by Acupuncturists Located above your head (above the crown)
Chakra Crown (CH Crn) (R)	Used by Acupuncturists Located same as BioTerminal Crown –Top of Head
Chakra 3rd Eye (CH 3rd) (R)	Used by Acupuncturists located same as BioTerminal Head -Forehead
Chakra Throat (CH Thr) (R)	Used by Acupuncturists Located same as BioTerminal® Neck -Thyroid
Chakra Heart (CH Ht) (R)	Used by Acupuncturists located same as BioTerminal® Chest- Heart
Chakra Solar Plexus (CH Solrp) (R)	Used by Acupuncturists located sames as BioTerminal® Abdomen
Chakra Sacrum (CH Sac) (R)	Used by Acupuncturists located same place BioTerminal® Pelvic
Chakra Root (CH Rt) (R)	Used by Acupuncturists located same place BioTerminal® Base
Blue Relax(R)	Scenar Russian technology Dr. Tennant wanted to add for those that use the SCENAR
Modulate	Scenar Russian technology Dr. Tennant wanted to add for those that use the SCENAR
Stimulate (Scenar FM/Var)	Scenar Russian technology Dr. Tennant wanted to add for those that use the SCENAR
Acute (Freq 121, Damp 2, Mod 3:1)	Scenar Russian technology Dr. Tennant wanted to add for those that use the SCENAR
Power 1	Toning face light
Power Plus	Toning of muscles
Power Advanced	Toning Body
Power Intense	Toning Abs Muscles, Arms, Thighs

Program	Feature
PG 2500	Highest Hertz
HGH	Human Growth Therapy Hormone
Repetitive Stimulate (R-Stim)	Used by Physical Therapists-
Fast T (R)	Used by Physical Therapists- Fast Twitch Muscle Fibers
Slow T (R)	Used by Physical Therapists- Slow Twitch Muscle Fibers
Lung (LU) (R)	Use to target organ for treatment
Pericardium (PC) (R)	Use to target organ for treatment
Heart (HT) (R)	Use to target organ for treatment
Small Intestine (SI) (R)	Use to target organ for treatment
Triple Burner (TE) (R)	Represents three organs at same time for treatment
Large Intestine (LI) (R)	Use to target organ for treatment
Spleen (SP) (R)	Use to target organ for treatment
Liver (LR) (R)	Use to target organ for treatment
Kidney (KI) (R)	Use to target organ for treatment
Bladder (BL) (R)	Use to target organ for treatment
Stomach (ST) (R)	Use to target organ for treatment
Advanced Mode (AVA) (where you program additional frequencies)	How to program a frequency, use this mode first before memorize mode
Memorized Mode (VAZ 4)	memorized mode
Memorized Mode (VAZ 3)	memorized mode
Memorized Mode (VAZ 2)	memorized mode
Memorized Mode (VAZ 1)	memorized mode
ER	Use in an Emergency, full Power, double select middle button